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Material:

Eldy Plus®

1. Kontaktdaten

 DentalPlus GmbH
Kohlgrub 5
83122 Samerberg
Tel.: 08032 9892007
Fax: 08032 9882790
E-Mail: info@dentalplus.info
Internet: www.dentalplus.info

2. Eigenschaften des Eldy Plus®

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Instruction manual injection molding

Material:

Eldy Plus®

1. Contact details

 DentalPlus GmbH
Kohlgrub 5
83122 Samerberg
Germany
tel.: +49 8032 9892007
fax: +49 8032 9882790
mail: info@dentalplus.info
url: www.dentalplus.info

2. Characteristics/properties of Eldy Plus®

- extreme break stability
- venerable
- high chemical bond to other plastics possible
- color fidelity and high translucency

3. Composition

Acrylic polymer based on methyl methacrylate (PETG)

4. Material-mechanical specifications

E- Modul	Ø 1.612 MPa	DIN EN ISO 20795-1:2013 compliant
Tensile strength	Ø 39,5 MPA	
Fracture stress	-	
Flexural strength	Ø 61,9 MPa	

5. Delivery forms of the cartridge

[See current injection molding catalog](#)

6. Indication


The material is suitable for producing

- Splints
- Partial and full dentures
- Crowns and bridges
- Framework for full veneers
- Model casting constructions
- Telescopic abutment constructions

Notes on Design

Minimum wall thickness:

- 0,6 mm (zervikal)
- 1,2 mm (okklusal)

8. Instructions for milling and processing

Eldy Plus® cartridges can be injected with our injection moulding device or with the injection moulding device from the company bredent.

To achieve an optimal injection moulding result, we recommend the use of our injection moulding equipment and the correct execution of the work instruction.

The restorations can be finished with conventional grinding and polishing agents and polished to a high gloss.

Processing instructions Eldy Plus®

- injection moulding temperature 270° C
- Pressure 7 - 8 bar
- Repressing time 60 seconds

(do not heat above 300 degrees, see safety data sheet)

You can find more detailed information at:

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9. Customizing / Repairs

Individual veneers should only be done as a full veneer. As a bonder we recommend visio.link® from bredent. Support in the cervical area and mechanical retention are beneficial. The veneering material is processed according to the manufacturer's instructions. Repairs can be carried out with cold polymers from AcrylX.

10. Cleaning

For gentle cleaning, we recommend conventional soapy water or detergent in an ultrasonic device or hand brush.

11. Attachment

For the attachment of restorations made out of **Eldy Plus®** all proven provisory and permanent cements are suitable.

12. Interactions

Interactions with other products commonly used in dentistry are not to be expected. If interactions cannot be ruled out in the mouth, the dental prosthesis must be removed

13. Storage

Eldy Plus® cartridges should be stored in a dry place at room temperature and protected from direct sunlight

14. Expiry date

See label on the packaging

15. First aid measures

Due to the processing of **Eldy Plus®** composite cartridges and its properties and in accordance with the instructions for use no health hazards are to be expected. If due to special circumstances any reactions take place, please note the following instructions:

- After inhalation: Seek medical advice if symptoms occur.
- After skin contact: Rinse with water and seek medical advice if symptoms occur.
- After eye contact: In case of mechanical irritation of the eyes, rinse with water and seek medical advice.
- After ingestion: If symptoms occur, seek medical advice.

16. Warnings

- **The use of denture finishers can lead to embrittlement (fracture) of the material if the corresponding ingredients are used.**
- Prolonged contact with high-proof alcohol solutions should be avoided. We recommend Bacoban for disinfecting plastics.
- Please note that improper use of milling tools and polishing brushes can cause the material to overheat. This can lead to damage to the material.